

For In Vitro Diagnostic Use
For Over-the-Counter (OTC) Use



Cue Health Monitoring System
Get Started Here




Cue Health




Set Up Your Cue[®] Account

Download the Cue Health App from the Apple App[®] Store[®] or Google Play Store[™] to your mobile smart device. Go to www.cuehealth.com for a list of compatible devices. Follow the on-screen instructions to create your account and pair your Cue Health Monitoring System (Reader) to your mobile smart device.



Connect Reader to power

When prompted on the Cue Health App, connect your Reader directly to power using the supplied Cue Power Adapter and USB Charging Cable. Only use the Power Adapter and USB Charging Cable supplied by Cue. Do not use if the Cue Reader is damaged.



Find Reader QR code

Rotate your Reader upside-down so the bottom of the Reader is facing up and the QR code is visible to you.

The QR code is located on the bottom of the Reader.



Pair Reader to mobile device

When prompted on the Cue Health App, adjust the camera so the QR code is visible on-screen and inside of the capture outline.






When the QR code is recognized, the Reader is paired securely to the mobile smart device.

Cue Reader is ready to use

To begin using the Reader to run tests, open the Cue Health App and follow the on-screen instructions.

Cleaning and Disinfecting

The Cue Reader should be cleaned and disinfected after each use. Wipe down with germicidal wipes containing 0.55% sodium hypochlorite as the active ingredient.

-  Do not spray any cleaning solution directly onto the Cue Reader or into the cartridge port.
-  Do not put any part of the Cue Reader under water or any other liquid.
-  Do not attempt to clean any internal parts.
-  If part of your device becomes wet, wipe off all moisture and allow sufficient time for drying before use.
-  Only germicidal wipes should be used. Other products have not been tested and should not be used. Follow manufacturer's instruction for handling and storage of wipes.

Welcome to Cue®

Welcome

This guide will help you get started using the Cue Health Monitoring System (Reader) and the Cue Health App.

Intended Use

The Cue Health Monitoring System (Reader) is a table top analyzer intended to be used with the Cue Health App and Cue Tests that have been cleared or approved with this instrument.

Get Started

See opposite page for instructions on downloading the Cue Health App, creating your Cue Account, and setting up your Cue Reader.

Questions

Dial (833) CUE-TEST · (833) 283-8378 · www.cuehealth.com

Cue and the Cue logo are registered trademarks of Cue Health Inc.

Apple and App Store are registered trademarks of Apple Inc.

Google Play is a trademark of Google LLC.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Cue is under license. Other trademarks and trade names are those of their respective owners.

About the Cue Health Monitoring System (Reader)

Cue Reader

Your Reader pairs with the mobile smart device, accepts Cue test-specific Cartridges, and runs tests.



About the Cue Health Monitoring System (Reader)



Your Cue Reader communicates with the Cue Health App using secure BLUETOOTH® wireless technology.


A Cue Reader can be used with multiple mobile smart devices, but only one smart device at a time can connect to a Reader to run a test and receive test results.



Charging the Reader

Charging with USB

Your Cue Reader can be quick-charged directly using a cable through the USB port located on the Reader.




Indicator Lights

Reader lights communicate status

The Cue Cartridge Reader features 5 LED indicator lights located on the top of the Cue Cartridge Reader.

The lights provide quick feedback on the status of the Cue Reader together with the Cue Health App.

The lights display unique patterns to communicate charging status, battery level, prompts for testing actions, test status, and test errors.



Battery Level

When the Reader is connected to power, the lights will sequentially illuminate then flash to indicate the battery charge level.



Lights illuminate sequentially



Battery is 80% charged



Battery is 20% charged

Prompts for an Action

While running a test or after a test has completed, the indicator lights will sweep from top to bottom or in reverse



Insert Cue Sample Wand with collected specimen into the Cue Test Cartridge



Remove Cue Test Cartridge (with Sample Wand inside) from the Cue Reader

Alerts

When connecting to the Cue Reader or at the start of a test, the indicator lights will flash in different patterns



Lights flash when connected



Lights flash when the Cue Test Cartridge or the Cue Sample Wand is inserted



Error. See the Cue Health App for details